


Name \_\_\_\_\_

### Agarose Gel Interpretation

FS 362


1. A PCR assay was designed to detect the *Salmonella* spp. gene *invA*. The *invA* PCR product is 375 basepair. The lanes (wells) are as follow: 1) Marker/standard, 2) *E. coli*, 3) *Listeria monocytogenes*, 4) *Bacillus subtilis* 5) Negative control 6) Salmonella, 7) unknown, 8) Marker



1. What are the purposes of lanes 1 & 8
2. What is the role for the reactions in lanes 6 & 7
3. Interpret lane 8
4. Identify 3 ways to improve this assay

2. The PCR assay below is a multiplex PCR designed to differentiate among STEC *E. coli*. "Shiga toxinogenic *Escherichia coli* (STEC) is an important cause of gastrointestinal disease in humans, particularly since such infections may result in life-threatening sequelae such as hemolytic-uremic syndrome (HUS). There is an increasing demand for improved diagnostic procedures for the detection of STEC in fecal samples and, in particular, in foods such as meat and dairy products. It has been recognized for a number of years that STEC strains causing human disease may belong to a very broad range of O serogroups. However, many of the STEC strains found in the gastrointestinal tracts of domestic animals (the principal source of human infections) may have a low degree of virulence in humans. These strains are less likely to produce putative accessory virulence factors such as intimin (encoded by *eaeA*) and the plasmid-encoded enterohemolysin (encoded by enterohemorrhagic *E. coli* (EHEC) *hlyA*). Within the human disease-associated strains, those producing Shiga toxin type 2 (Stx2, encoded by *stx<sub>2</sub>*) appear to be more commonly responsible for serious complications such as HUS than those producing only Shiga toxin type 1 (Stx1, encoded by *stx<sub>1</sub>*). We have developed multiplex PCR assays for the simultaneous detection of *stx<sub>1</sub>*, *stx<sub>2</sub>*, *eaeA*, and EHEC *hlyA*."

**Paton, A. W. et al. 1998. J. Clin. Microbiol. 36(2):598-602**



*stx<sub>1</sub>* = 534 bp

*stx<sub>2</sub>* = 384 bp

*eaeA* = 255 bp

EHEC *hlyA* = 180 bp

**Interpret each lane and state which isolates are most and least likely to cause disease.**

Name \_\_\_\_\_

**FS 362  
Primer Design Group Project**

Your goal is to design primers to amplify an approximate 600 base pair fragment of the *Salmonella* gene *invA*. *InvA* is essential for *Salmonella* spp. to enter epithelial cells. It is a member of a family of proteins involved in either flagellar biosynthesis or the secretion of virulence determinants by a number of plant and mammalian pathogens. This gene has been shown to be a suitable target for identification of *Salmonella* spp.

**The rules for Primer Design:**

- 1) Melting temperature should be between 58° and 60°C; melting T for forward and reverse primers should be the same
- 2) Avoid long “runs” of A’s, T’s, G’s, or C’s
- 3) Primers are always written 5’ to 3’, regardless whether it is the forward or reverse primer
- 4) Primers are between 18 and 24 bases long
- 5) Among the 5 bases on the 3’ end of the primer, 3/5 must be a G or a C.
- 6) Avoid palindromes which may result in hairpin loops
- 7) DNA Polymerase recognizes the 3’ free hydroxyl group.

Forward Primer:

5’ \_\_\_\_\_ 3’

Reverse Primer:

5’ \_\_\_\_\_ 3’

SOURCE: *Salmonella enterica* subsp. *enterica* serovar Saintpaul str. SARA29

LOCUS: *invA*

GENE: 2058 bp

MOLECULE: DNA

1 gtgctgcttt ctctacttaa cagtgcctgt ttacgacctg aattactgat tctggtacta  
61 atggatgatga tcatttctat gttcgtcatt ccattaccta cctatctggt tgatttctctg  
121 atcgactga atatcgtact ggcgatattg gtgtttatgg ggtcgttcta cattgacaga  
181 atcctcagtt tttcaacggt tcctgcggta ctgttaatta ccacgctctt tcgtctggca  
241 ttatcgatca gtaccagccg tcttatcttg attgaagccg atgccggtga aattatcgcc  
301 acgttcgggc aattcgttat tggcgatagc ctggcggggtg gttttgttgt cttctctatt  
361 gtcaccgtgg tccagtttat cgttattacc aaagggtcag aacgcgtcgc ggaagtcgcg  
421 gcccgatattt ctctggatgg tatgcccggg aaacagatga gtattgatgc cgatttgaag  
481 gccggtatta ttgatgcgga tgccgcgcgc gaacggcgaa gcgtactgga aagagaaagc  
541 cagctttacg gttcctttga cggcgatg aggtttatca aaggtagcgc tattgcccgc  
601 atcattatca tctttgtgaa ctttattggc ggtatttcgg tggggatgac ccgccatggt  
661 atggatttgt cctccgccct gtctacttat accatgctga ccattggtga tgggtctgtc  
721 gccagatcc ccgcattggt gattgagatt agtgccggtt ttatcgtgac tcgcgtaaat  
781 ggcgatagcg ataatatggg acggaatatc atgacgcagc tgttgaacaa cccatttgta  
841 ttggttgtta cggctatttt gaccatttca atgggaactc tgccgggatt cccgctgccg  
901 gtttttggtta ttttatcggg ggttttaagc gtactcttct attttaaatt ccgtgaagca  
961 aaacgtagtg ccgccaacc taaaaccagc aaaggcgagc agccgctcag tattgaggaa  
1021 aaagaagggt cgtcgttagg actgattggc gatctcgata aagtctctac agagaccgta  
1081 ccggtgatat tacttgtgcc gaagagccgg cgtgaagatc tggaaaaagc tcaacttgcg  
1141 gagcgtctac gtagtcagtt ctttattgat tatggcgtgc gcctgccgga agtattgtta  
1201 cgcgatggcg agggcctgga cgataacagc atcgtattgt tgattaatga gatccgtggt  
1261 gaacaattta cggctctattt tgatttgatg cgagtggtaa attattccga tgaagttggt  
1321 tcctttggca ttaatccaac aatccatcag caaggtagca gccagtattt ctgggtaacg  
1381 catgaagagg gggagaaact ccgggagctt ggctatgtgt tgcggaacgc gcttgatgag  
1441 ctttaccact gtctggcggg gacgctggcg cgcaacgtca atgaatattt cggatttcag


1501 gaaacaaaac atatgctgga ccaactggaa gcgaaatttc ctgatttact taaagaagtg  
1561 ctcagacatg ccacggtaca acgtatatct gaagttttgc agcgtttggt aagcgaacgt  
1621 gtttccgtgc gtaatatgaa attaattatg gaagcgctcg cattgtgggc gccaagagaa  
1681 aaagatgtca ttaaccttgt ggagcatatt cgtggggcaa tggcgcgтта tatttgccat  
1741 aaattcgcca atggcgggtga attacgagca gtaatggtat ctgctgaagt tgaggatggt  
1801 attcgcaaag ggatccgtca gacctctggc agtaccttcc tcagccttga cccggaagcc  
1861 tccgctaatt tgatggatct cattacactt aagttggatg atttattgat tgcacataaa  
1921 gatcttgtcc tccttacgtc tgtcgatgtc cgtcgattta ttaagaaaat gattgaaggt  
1981 cgttttccgg atctggaggt tttatctttc ggtgagatag cagatagcaa gtcagtgaat  
2041 gttataaaaa caatataa

Name :

### PFGE Gel Activity


You and your group are currently working for PulseNet Canada, investigating an outbreak related to Tim Horton's custard-filled donuts and *Staphylococcus aureus*. From January to June 2011, there were 19 reported illnesses with 5 hospitalizations and 1 mortalities. No illnesses were reported from July to September. From October to December 2011, an additional 8 illnesses were recorded.

You and your team have been given the agarose gels resultant of separate PFGE runs from each month. Below is a summary of the merged gels from the first occurrence (6 months):



Question 1: How many different strains were present during the first 6 months of the outbreak?


Following the interpretation of the first 6 months, you are provided with the data from the last three months (October to December).



Question 2: How many different strains were present during the last 3 months of the outbreak?

Question 3: Are there any similar strains from the first 6 months compared to the last 3 months?

Following the occurrences, samples were taken from all the Tim Horton's restaurants in the surrounding vicinity of living quarters for the reported cases and the victims.



Question 4: How many outbreaks are related to each store?

Store	Outbreak Number
Kamloops, B.C.	
Kelowna, B.C.	
Edmonton, A.B.	
Calgary, A.B.	
Calgary, A.B.	
Kimberly, B.C.	
Medicine Hat, A.B.	
Red Deer, A.B.	
Lethbridge, A.B.	



First 6 months

Last 3 months

Stores

